

Dawascope PIMS Installation Guidelines

Thank you for choosing Dawascope PIMS! We are truly delighted to have you as our valued partner and are committed to making your experience seamless. Dawascope PIMS is designed to be easy to install, taking less than 30 minutes when you follow the step-by-step instructions below. With this guide and our additional resources, you'll be set up in no time.

Need Help? We're Here for You!

If you encounter any challenges during or after installation, please don't hesitate to reach out to us at **+254712979721 (Call, WhatsApp)**. We've got you covered! Remember, unless you report an issue, we will assume all is well on your end.

Make sure to follow our WhatsApp Channel for communications, especially when we release software updates. <https://whatsapp.com/channel/0029Va18os75q08IJ80ped0D>

If you ever face any challenge, now or x years to the future, our capacity to help you depends on you having adhered to instructions provided. If you don't, in as much as we will sympathize with you, we will not be able to help you to the desired outcome.

Minimum Expectations

It is expected that you and your pharmacy team have the basic computer knowledge. How to register an account in any website, how to download anything from internet, how to create a folder in your own computer, how to move a file/folder from one location to another, what is right clicking, double-clicking... these are basic computer skills that by virtual of having ever interacted with a computer you should be having. We can only train you how to operate our software but not how to operate/use your computer.

Training Packages

To support you in getting the most out of Dawascope PIMS, we offer flexible training options tailored to your needs:

1. **Virtually assisted Self-Training** – Free
Access a wealth of pre-recorded video tutorials and a comprehensive user guide to train at your own pace. If you are not sure of the resource to refer to, you may call Dawascope team for assistance. Note that the assistance will involve being directed to the specific video to watch unless where it is a direct enquiry in need of a simple direct answer. Your cooperation is highly required. There is no limit to the direct answer queries you may make.
2. **Virtual Training** – Kes 2000/day
If you would rather have one-on-one comprehensive virtual training instead of using the available free training resources.
3. **On-Premises Installation** – Kes 500 to 3500 (plus return ticket from Nairobi).
This is where you want someone to come to your pharmacy to help with installation only (no training)
4. **On-Premises Training** – Kes 2000 to 3,500/day (plus return ticket from Nairobi)
A Dawascope expert will visit your location to provide hands-on installation and training.

Note: Charges for option 2 to 4 must be paid in full at the time of booking.

Virtually assisted installation (not training) is free, and it is your right to ask for it. But it may consume more of your time. And we expect you to have the basic computer skills for us to guide you effectively. But before contacting us, at least ensure you have read this entire guide. When we direct you to a specific section of the installation guide, we expect you to know where to go, because you will have read it.

Self-Training Resources

We've created detailed resources to guide you through installation and daily operations:

- **Video Tutorials:** Watch step-by-step guides here: [Dawascope YouTube Channel](#).
- **User Guide:** Available within the system or via your pharmacy Dawascope account at <https://dawascope.co.ke/help-content>.

Why Do We Charge for Training?

At Dawascope, we are committed to making our systems accessible and affordable for all users. To keep Dawascope PIMS subscription costs low and provide free trial, we have designed the system to be 100% self-service and provided extensive free resources such as pre-recorded training videos and helpful guides. These resources empower you to get started with minimal support.

However, for users who prefer personalized, hands-on assistance, we offer additional training services. The fees for these services help us cover the time and resources required to provide in-depth, one-on-one guidance while ensuring the system remains cost-effective for everyone.

We understand that every user's needs are different, and we are here to support you in the way that works best for you—whether through our free self-training materials or tailored training packages.

FOR SEAMLESS INSTALLATION, OPEN THIS GUIDE IN THE COMPUTER YOU ARE INTENDING TO INSTALL THE DAWASCOPE PIMS

NB: Your computer antivirus settings might block the downloads. Ensure you override or change them where necessary.

System Requirements

If your facility has more than one computer, choose one which will serve as the Server.

The computer should meet the following minimum requirements.

- a. Windows 8/10/11 64-bit operating system. If your installed OS is 32-bit, you can download a 64-bit operating system from Microsoft and install.
- b. Microsoft Visual C++ 2015-2022 Redistributable x64 [if your computer does not have this, download and install from https://aka.ms/vs/17/release/vc_redist.x64.exe
- c. WebView 2 Runtimes; if your computer does not have this, download and install from <https://msedge.sf.dl.delivery.mp.microsoft.com/filestreamingservice/files/cea4cf18-02a7-45d4-82ce-c9cbdd22da22/MicrosoftEdgeWebView2RuntimeInstallerX64.exe>
- d. Minimum hardware specifications: 4GB RAM and Core i3 processor (higher specs recommended).

Even if your computer does not meet the last requirement, you may still install though it may be slow in performance.

Ensure that you are not using VPN, and your PROXY setting is off.

Your computer should also be connected to a stable internet connection. If your system does not connect to Wi-Fi, you can use USB tethering

Downloading the Installation files

ARE YOU SURE YOU HAVE INSTALLED THE ABOVE SYSTEM REQUIREMENTS? DO NOT CONTINUE IF YOU HAVEN'T.

[Click this link to download the Dawascope installation files into a zip file.](#)

Once downloaded, copy the Dawascope zip file into a new folder

In the new folder, select the Dawascope zip file and extract it. This will create two folders, server and client folders.

NB: Dawascope PIMS is a dual system consisting of two software components: Dawascope_server and Dawascope_pims. Together, these two function seamlessly as a unified system.

Dawascope_server should only be installed on one computer. Dawascope_pims is the one to be installed on all your computers within your one pharmacy.

Ensure you finish the whole set up as described here in one sitting. Do not proceed past this point if you cannot complete the installation and setup. You must maintain internet connection throughout.

Installing Dawascope Server

Choose the computer that will act as the server for your other pharmacy computers and implement this installation on it.

Ensure the system requirement b and c listed above are installed

Ensure the date and time of your computer is correct, and ensure it remains correct. The clock should be showing the current time, the time zone should be UTC + 03:00 Nairobi and the region should be Kenya and region format English Kenya.

Open the server folder you have extracted above

Right click on the Dawascope_server_setup application file and select Run as Administrator Install in the suggested folder and add a desktop shortcut.

After the installation, launch the software at the finish.

A black application terminal window will open, give it time to set up all the necessary files. **Make sure you do not click on this terminal window. If you do, you will pause the process, and you may not know.**

After all files are set, the server will start running and it will open a second window indicating **Your Server is Up and Running, kindly open the twin application Dawascope PIMS to utilize the service.**

Do not close these two Server windows. You can minimize them.

Installing Dawascope_pims

Open the client folder you have extracted above

Right click on the Dawascope_pims_setup application file and select Run as Administrator Install in the suggested folder and add a desktop shortcut.

After the installation, launch the software at the finish.

The Dawascope server should be running.

Registering Your Dawascope software

Launching Dawascope_pims will open the software registration page. (if it opens the sign in page before activation, close it and open using the desktop shortcut you have created)

Use the unique key you received in your email to activate your software. This key can only be used to activate one software. Therefore, ensure that the computer you are working on will permanently be your Server, and it will always be at the pharmacy.

Input the key and submit and your software is activated.

Then set up the invoice number prefixes for the modules you have subscribed to. The system is set up to generate a unique number for every sale you will make. You are required to set the prefix for this number that will uniquely identify it with your organization. For instance, when you sell, the software will generate a receipt number in numerals 457398. Your prefix will be attached to it to identify it with your pharmacy. XYZ457398.

You will then set the percentage profit margin. The system would set the selling price using this profit margin if you chose to move products automatically. But you can set different selling prices or update the set prices anytime you wish without limits. If you enter 5%, it means the selling price of an item bought at Kes 100 will be 105. If you set 30, the selling price will be 130. If you don't understand this, but only about calculating a markup of 1.33, 1.35.... 1.33 means 33% profit margin.

Then you will set the expiry timeline. This is the number of days that, if the expiry date of an item falls within the set days, that item will be marked as short expiry and be notified to you.

Then you will update your catalog which will take less than 5 minutes. Ensure your internet connection is reliable.

You are all set to start using your Dawascope PIMS system.

Multiple Computers

If you have more than one computer in your pharmacy, copy the Dawascope_pims file to the other computer(s) and install it there. Do not install the Dawascope_server on more than one computer.

Ensure all the computers in your pharmacy are connected to the same network (internet connectivity not required)

First Sign In

Sign in using the credentials **username: admin** and **password: admin**. You are advised to change the password. Go to Settings then Account settings to change the password. This is the super user and therefore ensure you secure your credentials. DO NOT FORGET THE ADMIN PASSWORD. IF YOU DO, IT WILL REQUIRE MANUAL INTERVENTION FROM DAWASCOPE TEAM. FACILITATION CHARGES WILL APPLY

Help page

Once logged in, click on the [?] icon on the top bar left to [Access Pharmacy] button to access the help page for more details on using your software. You can also access the help page from your online portal using the phone for easier reference. You can also watch YouTube videos.

Creating staff accounts

You are advised to restrict use of the admin account for administrative purposes only.

Even if you are the only person in the pharmacy, do not use the admin account to transact. Create your staff account. Follow the guidelines on the help page on how to create and manage staff (employee's) account.

We also notify you via help page when we release updates (you will need to sync catalog whenever you need to see the updated notification)

Subsequent use of your Dawascope PIMS

Whenever you want to use your software, **Right Click** on the Dawascope Server shortcut on the desktop and select Run as Administrator and accept the prompt. You should have created the shortcut during installation as guided above. **If you do not run the server with Administrator privileges, some components will not work.**

Avoid clicking on the server terminal window.

Then double click on the Dawascope PIMS shortcut on the desktop or via the windows program list. The PIMS will open and connect to the Server. It will display the login page.

Closing Dawascope Server

You should follow the correct procedure when shutting down the Dawascope Server. Do not use the normal X button at the top right corner of the running window. **We cannot stress this enough.**

Open the help page and you will see the guide on how to shut down the Server on the first page.

[CTRL+C]

Backup of your Database

This is not an automated service; it must be triggered manually. When you create your staff account, give them the permission to do the database backups. You are advised to do backup daily when you are closing your business. If you operate 24 hours, decide when to do the backups. You can do as often as you deem fit, but once a day is sufficient. Do not attempt to interfere with the backup files. If your system experiences any problems, we will use the most recent and valid backup

file to roll you back to that date. If you don't provide us with proper backup file, we will not be able to assist you in recovering your system.

When the backup is triggered, a protected file is created in your downloads folder with the name database_backup_1234567890. Do not open this file. If you alter it, you will invalidate it.

Note that this backup is within your computer. If your computer is stolen, you will lose it. If you use your computer for other purposes that predispose it to malware attacks, and it is attacked with malware, the backup may also be corrupted. You are therefore advised to occasionally copy the backup file to external private storage and/or send to your secure email address.

Follow these simple steps to configure backup. Execute this only in the computer running the Dawascope Server.

- i. Click this link to download Postgres application
<https://get.enterprisedb.com/postgresql/postgresql-15.12-1-windows-x64.exe>
- ii. Once downloaded, RIGHT CLICK on the downloaded file and select Run as Administrator.
- iii. Follow the prompts to install the PostgreSQL, do not change the installation directory, in the next page, ensure all components are selected, do not change data directory,
- iv. Then you will be prompted to enter password. Input XYZ (all uppercase) as the password.
- v. Do not change the port, leave it 5432
- vi. Locale should remain [Default locale]
- vii. You will then see installation summary, click next
- viii. Ready to Install, click next, and it will start installing.
- ix. After it has installed, it will show the final window indicating that it has finished installation. Before clicking finish, **UNMARK the Stack Builder.....** Then click finish
- x. After you have done the above, click this link to download the Postgres environment variable setter <http://dt.dawascope.com/DawascopePSQLEnVariable.bat>
- xi. Once downloaded, RIGHT CLICK on the DawascopePSQLEnVariable file and select Run as Administrator. It should take 1-5 seconds to complete the task
- xii. Restart the computer, rerun the Dawascope system (You must Run as Administrator) and trigger the Database Backup. You can see Download Backup button at Settings>General Settings when logged in as admin. Staff can use the Backup Local DB under Cashier menu. You should see a file generated.

Printer setup

Only 80 mm Thermal Receipt Printer with the following specifications is supported.

1. Interfaces: USB
2. Command: ESC/POS

When buying such a printer, ensure it has the above on its specification tag. Avoid Epson printer due to complicated settings involved.

Go to the help page, click on Accessory set-up for the guide on how to set up your printer.

Always remember that this system is made up of two applications and both applications must be running.

Protect your computer from malware. Avoid strange flash-disks, downloading strange files. Limit the use of this computer to Dawascope System only.

Your Success is Our Priority

At Dawascope, your success is our success. Whether you're installing the software, training, or operating day-to-day, we're here to assist. Reach out anytime for support or training arrangements.

Thank you once again for choosing Dawascope PIMS. We're thrilled to have you on board and look forward to empowering your pharmacy operations!

Warm regards,

The Dawascope Team